

## BAG – Batch Monitoring Systems

**STEAM**



### **BAG-RePac 121 °C / 134 °C / 18 Min.**

Batch Monitoring Systems for monitoring steam sterilization cycles in pre-vacuum sterilizers at 121 °C or 134 °C. Chemical indicator (ISO 11140-1: Class 6) integrates the parameters steam, temperature and time. Detection of residual air and leakage problems. To be used with RePac test device.

**Packaging unit:** 300 tests/box

<b>BAG-RePac 121 °C</b>	<b>Order-No. 7541</b>
<b>BAG-RePac 134 °C</b>	<b>Order-No. 7542</b>
<b>BAG-RePac 18 (134 °C, 18 min.)</b>	<b>Order-No. 7553</b>
<b>BAG-RePac test device</b>	<b>Order-No. 7540</b>

**STEAM**



### **BAG HeliPac 121 °C / 134 °C / 18 Min.**

Batch Monitoring Systems for monitoring steam sterilization cycles in pre-vacuum sterilizers at 121 °C or 134 °C. Chemical indicator (ISO 11140-1: Class 6) integrates the parameters steam, temperature and time. Detection of residual air and leakage problems. Further development of BAG-RePac with helix type test device.

**Packaging unit:** 1,200 tests/box incl. helix

<b>BAG-HeliPac 121 °C</b>	<b>Order-No. 7571</b>
<b>BAG-HeliPac 134 °C</b>	<b>Order-No. 7572</b>
<b>BAG-HeliPac 18 (134 °C, 18 min.)</b>	<b>Order-No. 7573</b>

Contents of the system:

- HeliPac test device (helix)
- HeliPac indicators, which ensures secure results with the approved purple to green colour change
- LEAK test indicators for self-testing and security of the system

-----  
Also available: Packaging units with 600 tests/box incl. helix  
-----

**Refill package:** 300 tests/box without helix

<b>BAG-HeliPac 121 °C</b>	<b>Order-No. 75711</b>
<b>BAG-HeliPac 134 °C</b>	<b>Order-No. 75721</b>
<b>BAG-HeliPac 18 (134 °C, 18 min.)</b>	<b>Order-No. 75731</b>

**STEAM**

## BAG-Helix T

For table-top sterilizers only.

Single-use Batch Monitoring Systems for monitoring steam sterilization cycles in pre-vacuum sterilizers at 121 °C or 134 °C (e.g. type B acc. To EN 13060). Chemical indicator (ISO 11140-1: Class 5) integrates the parameters steam, temperature and time. Detection of residual air and leakage problems.

Product features:

- Helix test device acc. to EN 867-5
- Chemical indicators which ensures secure results with the approved purple to green colour change
- Safety: To be opened only after sterilization
- Single-use product, no preassembling necessary

**Packaging unit:** 50 tests/box

**BAG-Helix T**

**Order-No. 7576**

02/09